



# A new class of compact ultrasound

Philips CX50 CompactXtreme ultrasound  
system specifications

**PHILIPS**  
sense and simplicity

# Table of contents

<b>1. Introduction</b>	<b>3</b>	<b>Linear array</b>	<b>11</b>
Key advantages	3	L12-3 broadband linear array	11
1.1 Applications	3	<b>Sector arrays</b>	<b>11</b>
1.2 Optional portability	3	S5-1 sector array with PureWave technology	11
<b>2. System overview</b>	<b>4</b>	S8-3 broadband sector array	11
2.1 System architecture	4	S12-4 broadband sector array	11
2.2 Imaging modes	4	X7-2t xMATRIX array with PureWave technology	11
2D grayscale	5	<b>Non-imaging</b>	<b>11</b>
M-mode	5	D2cwc CW transducer (Pedoff)	11
Anatomical M-mode	5	D5cwc CW transducer (Pedoff)	11
Tissue Doppler imaging	5	5.3 Transducer application guide	12
Pulsed wave Doppler	5	<b>6. Measurements and analysis</b>	<b>13</b>
Continuous wave Doppler	5	6.1 Comprehensive measurement tools	13
Spectral Doppler	5	6.2 QLAB advanced quantification	13
Color Doppler	6	Region of Interest (ROI) quantification plug-in	13
Contrast imaging	6	Tissue Motion Quantification (TMQ) plug-in	13
Tissue harmonic imaging	6	Strain Quantification (SQ) plug-in	13
Color Power Angio imaging	6	Intima Media Thickness (IMT) measurement plug-in	13
Freehand 3D	6	MicroVascular Imaging (MVI) plug-in	13
<b>3. System controls</b>	<b>7</b>	GI 3D Quantification (GI 3DQ) plug-in	13
3.1 Advanced imaging controls	7	<b>7. Physical specifications</b>	<b>14</b>
3.2 Control panel and user interface	7	System dimensions	14
<b>4. Workflow</b>	<b>8</b>	Physical features	14
4.1 Display annotation	8	Mobility cart	14
4.2 Image presentation	8	Travel case	14
4.3 Cineloop review	8	Power requirements	15
4.4 Connectivity	9	Power management	15
4.5 Ergonomics	9	Environmental	15
4.6 Stress echo protocols	9	Electrical safety standards	15
4.7 Integrated intervention	9	Safety requirements	15
<b>5. Transducers</b>	<b>10</b>	Maintenance and serviceability	15
5.1 Transducer selection	10		
5.2 PureWave transducer technology	10		
<b>Curved arrays</b>	<b>10</b>		
C5-1 curved array with PureWave technology	10		
C8-5 broadband curved array	10		
C9-3v curved array with PureWave technology	11		

# 1. Introduction

The new CX50 system is built on a new platform and architecture in a small, compact design that is ideal for taking premium performance anywhere you need it. PureWave, a clinically proven imaging technology previously available only on Philips premium systems, captures a broad band of tissue information for exceptional clarity and information. Tissue Specific Imaging presets offer a new level of exam automation – transducers are optimized by exam type, providing excellent images with minimal adjustment. Highly configurable portability enables you to study the most difficult-to-image patients across a variety of clinical needs.

The CX50 is fully configurable, allowing you to select imaging capabilities, transducers, and clinical analysis to support your exam needs. Add supported capabilities at any time as well as upgrades when they become available.

## Key advantages

- Premium performance anywhere you need it
- Exceptional results on technically difficult patients
- Support for cardiology, general imaging, vascular, and Ob/Gyn exams

## 1.1 Applications

Adult cardiac
Adult transesophageal
Stress echo
Abdominal
Pediatric
Vascular
Obstetrical
Gynecological and fertility
Small parts
Breast
Musculoskeletal
Emergency medicine
Regional anesthesia
Intervention
Contrast

## 1.2 Optional portability

Specially designed cart
Wheeled travel case
Additional AC adapter



## 2. System overview

Philips proprietary technologies are an integral part of the CompactXtreme platform, and provide the basis for its extensive range of imaging capabilities.

### 2.1 System architecture

- Next generation all-digital compact broadband beamformer with pulse shaping capability
- High resolution A/D conversion with 170 dB full time system dynamic range
- 20,000 digitally-processed channels
- Supports PureWave technology and high-density array configurations
- Multi-variate harmonic imaging including pulse inversion processing
- One-touch 2D optimization with broadband frequency compounding
- SonoCT real-time beam-steered compound imaging
- Advanced XRES adaptive image processing
- Continuously variable steering in 2D, color Doppler, and Doppler modes
- iSCAN one-touch intelligent optimization for 2D and Doppler
- Active native data manipulation
- Tissue specific imaging presets
- Gray shades: 256 levels (8 bit) in 2D, M-mode, and Doppler

### 2.2 Imaging modes

- 2D
- M-mode
- Anatomical M-mode
- Color M-mode
- Color Power Angio (CPA) imaging
- Directional CPA
- Pulsed wave (PW) Doppler
- HPRF PW Doppler
- Continuous wave (CW) Doppler
- Freehand 3D imaging
- QLAB advanced quantification software
- Invert and color invert
- Color compare mode
- Dual mode
- Duplex for simultaneous 2D and Doppler
- Triplex for simultaneous 2D, Doppler, and color or Color Power Angio



The CX50 system combines premium imaging technologies and advanced optimization for exceptional image clarity results.

- 2D and flow optimization signal processing
- Intelligent Doppler – automatically maintains pre-selected 0/60 degree flow angle
- Live compare
- Tissue harmonic imaging (THI)
- Reconstructed zoom with pan (read zoom)
- High definition write zoom
- Trapezoidal imaging
- Pulse inversion harmonic imaging
- Contrast imaging
  - Cardiac left ventricular opacification (LVO)
  - General
  - Superficial
- Adaptive Doppler
- Adaptive color Doppler
- Color tissue Doppler imaging
- Pulsed wave tissue Doppler imaging
- Active native data (allows manipulation of raw image data)
- SmartExam system-guided protocol capability

### 2D grayscale

- Smart TGC: pre-defined TGC curves optimized for consistently excellent imaging with minimal TGC adjustment
- User adjustable LGC control
- High definition zoom concentrates all image processing power into a user-defined area of interest; possible to combine high definition zoom with pan zoom
- Cineloop image review
- Selectable 2D compression settings
- Sector size and steering control for sector and curved array image formats
- Dual imaging with either independent cineloop buffers or split screen imaging
- Live compare
- Chroma imaging with multiple maps
- 256 (8 bits) discrete gray levels
- 2D acquisition frame rate up to 755 frames/sec (dependent on field-of-view, depth and angle)

### M-mode

- Available on cardiac imaging transducers
- Selectable sweeping rates
- Time markers: 0.2 and 1.0 seconds
- Acquisition zoom capability
- Selectable display format prospective or retrospective (1/3-2/3, 2/3-1/3, side by side, full screen)
- Chroma colorization with multiple color maps
- Cineloop review for retrospective analysis

### Anatomical M-mode

- Available for cardiac imaging transducers
- Uses 2D image as a basis for M-mode analysis at a defined line, independent of transducer orientation
- Makes it easier to keep the M-mode line perpendicular to the anatomy, even in abnormally shaped or positioned hearts
- Provides data on direction, position and timing of any single echo received from any point of the tissue for M-mode analysis in any direction, for examining cardiac chamber diameters, LV regional wall motion, and location of accessory pathways
- Anatomical M-mode trace can be generated or modified post Freeze

### Tissue Doppler imaging

- Available on S5-1 transducer
- Allows high frame rate acquisition of tissue motion
- Color gain, TGC, and LGC
- 8 color maps
- Velocity (cm/s)

### Pulsed wave Doppler

- Available on all imaging transducers
- Adjustable sample volume size: 0.8-24.6 mm (transducer dependent)
- Simultaneous or duplex mode of operation
- Simultaneous 2D, color Doppler or CPA, pulsed Doppler
- iSCAN optimization automatically adjusts scale, baseline and Doppler gain

### Continuous wave Doppler

- Available on cardiac sector array transducers and non-imaging transducers
- Steerable through 90° sector
- Maximum velocity range: 20 m/sec (transducer dependent)

### Spectral Doppler

- Display annotations including Doppler mode, scale (cm/sec) Nyquist limit, wall filter setting, gain, acoustic output status, sample volume size, normal/inverted, angle correction, grayscale curve
- Angle correction with automatic velocity scale adjustment
- Adjustable velocity display ranges
- Normal/invert display around horizontal zero line
- Selectable sweep speeds
- Selectable low-frequency signal filtering with adjustable wall filter settings
- Selectable grayscale curve for optimal display
- Selectable Chroma colorization maps
- Selectable display format prospective or retrospective (1/3-2/3, 2/3-1/3, side by side, full screen)
- Doppler review for retrospective analysis of Doppler data
- 256 (8 bits) discrete gray levels
- Post-processing in PW frozen mode includes map, baseline, sweep speed, invert, compress, reject, and Chroma



### Color Doppler

- Adaptive mode adjusts Doppler frequency and sensitivity based on color ROI placement within image available on all imaging transducers
- Cineloop review with full playback control
- Advanced motion suppression with intelligent algorithms; adapts to various application types to selectively eliminate virtually all color motion artifact
- 256 color bins
- Trackball-controlled color region of interest: size and position
- Maps, filters, color sensitivity, line density, smoothing, echo write priority, color persistence, gain and baseline optimized automatically by exam type or is user selectable
- Velocity and variance displays
- Color invert in live and frozen imaging
- User selectable smoothing control
- User selectable persistence control
- Color/2D line density control

### Contrast imaging

- System optimized for detecting contrast agent signatures
- Contrast modes available on the S5-1, C5-1 and L12-3 transducers
- Pulse inversion contrast imaging available with XRES technology
- Power modulation (PM), pulse inversion (PI), coded harmonic, and flash contrast imaging modes
- Low MI mode
- Display timer
- Low MI color flow contrast
- S5-1 left ventricular opacification (LVO) for adult cardiology applications
- ECG/timed triggering
- QLAB ROI display

### Tissue harmonic imaging

- Second harmonic processing to reduce artifacts and improve image clarity
- Multivariate pulsing including patented pulse inversion phase cancellation technology for maximum detail resolution during harmonic imaging
- Available on all imaging transducers
- Extends high performance imaging capabilities to all patient body types
- Supports SonoCT (harmonic SonoCT) and XRES modes

### Color Power Angio imaging

- Highly sensitive mode for small vessel visualization
- Available on all imaging transducers
- Cineloop review
- Multiple color maps
- Individual controls for gain, PRF, baseline, filters, sensitivity, echo write priority and color invert
- Dynamic motion differentiation
- Adjustable CPA region of interest: size and position
- User-selectable persistence
- User-selectable blending
- Directional Color Power Angio (DCPA) mode

### Freehand 3D

- Qualitative grayscale volume acquisition supported on all imaging transducers
- Volume display with surface rendering (transparency, threshold, smoothing, brightness, and opacity controls)
- Multiplanar view display
- Specialized algorithms and maps maximize three-dimensional display
- Trim tools on both volume and multiplanar reconstructed (MPR) views
- Supported by SonoCT and XRES modes to reduce noise artifacts
- Resize control adjusts for different sweep speeds
- Sculpt/erase control of volume
- Advanced QLAB volume analysis tools:
  - 3DQA: MPR, thick slice, and advanced volume rendering capabilities
  - iSlice precision tomographic volume imaging capability

One-button controls are logically placed on the CX50 control panel for easy optimization during every exam.



# 3. System controls

The CX50 has sophisticated system controls to help you acquire the best possible data on your patients, including many one-button optimization controls that adjust thousands of system parameters.

## 3.1 Advanced imaging controls

### iSCAN image optimization

- One-touch image optimization
- In 2D mode, one button automatic adjustment of:
  - TGC and receiver gain to achieve optimal uniformity and brightness of tissues
- In Doppler mode, one button automatic adjustment of:
  - Doppler PRF based on detected velocity
  - Doppler baseline based on detected flow direction
  - Gain to achieve optimal brightness of spectral waveform
- Available on all imaging transducers
- Operates in conjunction with SonoCT and XRES imaging

### SonoCT real-time compound imaging

- Available on all curved transducers and linear array
- Eliminates virtually all clutter and artifact
- Automatic selection of the number of steering angles (up to 9) based on the user-selected resolution/frame rate (Res/Speed) condition
- Operates in conjunction with tissue harmonic imaging, volume modes, imaging, and duplex Doppler
- Operates in conjunction with XRES imaging
- Available in contrast modes

### Advanced XRES adaptive image processing

- Available on all imaging transducers
- Eliminates virtually all speckle noise and enhances border definition
- Available in contrast modes

### Expanded field of view

- Trapezoidal imaging
  - Expands field of view on linear array transducers up to 15° on each side in vascular and general imaging applications

### Active native data

- 2D image controls that can be changed in review include: gain (overall gain, TGC, LGC), compress, gray map, Chroma map, orientation (L/R, U/D), display zoom/pan, XRES
- PW and CW Doppler controls that can be changed include: gain, baseline, invert, angle correct, angle 60/0/60, sweep speed, grayscale and Chroma maps compress and reject, PW trace (High Q controls), display format
- Color image controls that can be changed in review include: gain, baseline, color map, invert, write priority, smoothing, suppress, variance, directional Color Power Angio
- Physio controls that can be changed: sweep speed, position, gain
- Can be acquired in prospective and retrospective timing sequences
- Images are acquired at acoustic data frame rate
- Available in cine loop and quick review modes

### Live compare

- Allows recall of current or previous exam image data for direct side-by-side comparison with current image data

## 3.2 Control panel and user interface

- Easy-to-learn graphical user interface
- Ergo-centric design of primary controls readily accessible and logically grouped
- Tri-state control panel lighting (active, available, and unavailable)
- Automatic ambient lighting sensing for optimal image viewing in both light and dark environments
- Dual function mode switch and independent gain controls for 2D, CPA, M-mode and color, PW, CW Doppler
- Eight-slide pot control adjustment of TGC curve
- Two-slide pot control adjustment of LGC curve
- iSCAN control for 2D/Doppler/color Doppler automatic optimization
- High definition/pan zoom control
- Freeze control
- Programmable print control
- Transducer selection and tissue specific imaging control
- Report and review controls
- Protocol selection control

# 4. Workflow

The CX50 system adapts to your workflow, whether you're in a critical care unit, at the bedside, in the emergency department, OR arena, or at a remote location. With easy-to-use tools designed for your needs, you're ready to scan wherever your patients are located.

## 4.1 Display annotation

- On-screen annotation of all pertinent imaging parameters for complete documentation, including transducer type and frequency, active clinical options and optimized presets, display depth, TGC curve, grayscale, color map, frame rate, compression map value, color gain, color image mode, hospital name, and patient demographic data
- User selectable display of patient birth date or user ID, institution name, and performed by
- Annotation data and patient name can be turned off (hidden) for generating images used in publication and presentation
- Scan plane orientation marker
- User selectable depth scale display
- Real-time display of mechanical index (MI)
- Real-time display of thermal index (TIb, TIc, TIs)
- Multiple trackball-driven annotation arrows
- Pre-defined body markers, supported in single and dual imaging formats
- Doppler baseline invert in live and frozen imaging
- TGC curve (user selectable On/Off display)
- TGC values (On/Off display)
- Tool Tips provides a brief description of the abbreviated on-screen image parameters
- Informative trackball arbitration prompts
- Thumbnail display of images printed/stored
- Calculations results and analysis labels
- User friendly menus that allow navigation to other analysis features
- Network and connectivity icons to allow instant feedback about network and printer conditions
- Cineloop frame display
- Cineloop bar with trim markers
- Prompt region for informational message display
- Protocol procedure list with status

## SmartExam protocols

- Exam guide with on-screen display
- Required views based on exam type
- Fully customizable protocol capability for clinical applications supported on the system with flexibility to conduct the examination protocol in any sequence
- Preset protocols for abdominal, vascular, and gynecological exams based on industry and accreditation guidelines
- Automatic launching of annotation and body marker icon on required views
- Automatic launching of calculations
- Ability to pause and resume SmartExam function at any time
- System analysis capabilities supported in all defined protocols
- Custom protocol transfer between CX50 systems

SmartExam's on-screen display provides immediate visibility of exam status.



## 4.2 Image presentation

- Up/down
- Left/right
- Multiple duplex image formats (1/3-2/3, 2/3-1/3, 50/50, and full screen)
- Depth from 1 cm to 30 cm (transducer dependent)

## 4.3 Cineloop review

- Acquisition, storage, and display in real time and duplex modes of up to three minutes in quick review of 2D and color images
- Dual imaging (single and dual buffer)



SmartExam protocols allow you to focus on your patients while confident you are capturing the required views.

#### 4.4 Connectivity

- Two USB ports on control panel
- 80 GB hard drive space
- Internal slot-load CD/DVD RW drive
- DICOM print, store, and storage commitment
- DICOM structured reporting for cardiac and obstetrics
- Performed procedure step (PPS)
- Modality worklist
- DICOM reader saved onto media
- Export data as PC-compatible or DICOM files
- Ethernet at 100 Mb/second
- Wireless “B and G” networking
- USB to serial converter adapter
- Support for optional small B/W and color printers

#### 4.5 Ergonomics

- Philips common user experience control panel with central track ball and easy-access mode keys
- Tri-state lighting allows immediate feedback of active and available controls in all modes
- High-resolution LCD display with wide viewing angle and automatic ambient light compensation
- Quick keys and active mode
- System-guided exam protocol capability

#### 4.3 Stress echo protocols

- Acquisition of single-frame or full-motion digital clips in any mode (including 2D, color Doppler, color TDI); type of image to be acquired may be changed on the fly by the operator as needed through pause protocol feature
- Gain Save adjusts automatically to different views
- Automatically saves your preferred control settings – such as MI (Mechanical Index), gain and depth – for each view while acquiring resting images
- At immediate post-exercise, system automatically retrieves saved settings for each view
- Allows different gain profiles for parasternal LAX and SAX views, AP4 and AP2 views
- Systole or full heart cycle acquisition
- Default stress protocols
  - May not be edited but may be used as the basis of a user-defined protocol
  - Factory-provided protocols include:
    - Two-stage exercise stress
    - Four-stage pharmacological stress
    - Three-stage exercise stress



- User-defined stress protocols
  - Utility for creation of user-defined protocols and editing of existing protocols for acquisition of stress and routine images, allows protocols to be defined to do any or all of the following:
    - Support between 1 and 8 stages
    - Support user-defined stage names
    - Support between 1 and 8 views per stage
    - Support user-defined view names
    - Prompt for a particular stage and view
    - Assign stage and view names
    - Set the number of cycles/beats for each view
    - Define prospective or multi-cycle/full disclosure acquisition
    - Save user-defined protocols within a preset
    - Save user-defined protocols to removable media for import onto other CX50 systems at the same software level
    - Modify protocols during use
    - Add stages at any point after the current stage
    - Change the name of a stage at any point up to acquisition of the first image of the stage
    - Save the modified protocol (it will not be automatically saved)

#### 4.7 Integrated intervention

- Integration capability with Philips Allura X-ray systems for controls and display
- Support for radiology, vascular, and cardiology interventional procedures

# 5. Transducers

The CX50 ultrasound system offers a full complement of transducer options that extend capabilities to meet a wide range of imaging needs.

## 5.1 Transducer selection

- Automatic parameter optimization of each transducer for exam type through Tissue Specific Imaging (TSI) software
- User-customizable imaging presets for each transducer
- Dedicated connector for continuous wave Doppler (pedoff) transducers
- Continuous dynamic receive focusing on all imaging transducers

## 5.2 PureWave transducer technology

- Proprietary combination of PureWave crystal, impedance matching layers, backing materials, micro-electronics, and cable design
- Achieves breakthrough acoustic broadband response and twice the efficiency of conventional technology for superb image quality and Doppler performance
- Designed for multi-variate harmonic imaging capabilities including pulse inversion and coded pulse sequencing
- Supports highly sensitive contrast enhanced ultrasound capabilities
- Available in sector, curved, tightly curved, and xMATRIX array configurations up to 12 MHz in frequency

## Curved arrays

### C5-1 curved array with PureWave technology

- 5 to 12 MHz extended operating frequency range
- High density curved array, 160 elements
- 2D, steerable PW Doppler, High PRF and color Doppler; and Color Power Angio, SonoCT, advanced XRES, and multi-variate harmonic imaging
- General purpose abdominal, obstetrical, gynecological, interventional, acute care, and regional anesthesia applications
- Contrast application
- Supports reusable, 4-angle, plastic biopsy guide (14-23 gauge)

### C8-5 broadband curved array

- 8 to 5 MHz extended operating frequency range
- High density curved array, 128 elements
- Field of view: 90 degrees
- Steerable pulsed Doppler, color Doppler, Color Power Angio, SonoCT, and advanced XRES imaging
- General purpose pediatric and vascular imaging
- Supports reusable, single angle, plastic biopsy guide (8.5FR and 14-23 gauge needles, 19 gauge not available)



CX50 supports a broad range of exam requirements with sector, curved, tightly curved, and linear array transducers

### C9-3v curved array with PureWave technology

- 9 to 3 MHz extended operating frequency range
- High density curved array, 128 elements
- End-fire sector, 11.5 mm radius of curvature, 130° field-of-view
- Steerable pulsed wave, High PRF, and color Doppler, Color Power Angio, SonoCT, advanced XRES, and multi-variate harmonic imaging
- Endovaginal applications
- Supports disposable, plastic biopsy guide (18 gauge) and stainless steel biopsy guide (16-18 gauge)

### Linear array

#### L12-3 broadband linear array

- 12 to 3 MHz extended operating frequency range
- Fine pitch, high resolution linear array, 160 elements
- 2D, steerable PW, High PRF, and color Doppler, Color Power Angio, SonoCT, advanced XRES, and multi-variate harmonic imaging
- Vascular, small parts, musculoskeletal, regional anesthesia, and acute care applications
- Contrast application
- Supports reusable, 3-angle, biopsy guide (14-23 gauge)

### Sector arrays

#### S5-1 sector array with PureWave technology

- 5 to 1 MHz extended operating frequency range
- Sector array, 80 elements
- 2D, steerable PW Doppler, CW Doppler, High PRF Doppler, color Doppler, tissue Doppler, advanced XRES and multi-variate harmonic imaging including contrast LVO
- Adult cardiac, general purpose abdominal, transcranial Doppler (TCD), acute care, and adult LVO applications
- Supports reusable, 3-angle, plastic biopsy guide (14-23 gauge) and 3-angle stainless steel biopsy guide (14-23 gauge)

#### S8-3 broadband sector array

- 8 to 3 MHz extended operating frequency range
- Phased array, 96 elements
- Scanplane aperture: 15.4; Field of view: 90°
- Steerable PW Doppler, CW Doppler, high-PRF Doppler, color Doppler, Tissue Doppler Imaging, harmonic imaging, and XRES imaging
- Pediatric and small adult cardiology and pediatric abdominal applications



The C5-1 transducer's design combines PureWave technology and ergonomics for outstanding image quality and user comfort, even on technically difficult patients.

#### S12-4 broadband sector array

- 12 to 4 MHz extended operating frequency range
- Phased sector array, 96 elements
- Scanplane aperture: 9.8; field of view: 90°
- Steerable PW Doppler, CW Doppler, High-PRF Doppler, color Doppler, Color Power Angio, XRES, and Tissue Doppler Imaging
- Pediatric and neonatal cardiology, pediatric abdominal, neonatal head applications

#### X7-2t xMATRIX array with PureWave technology

- 7 to 2 MHz extended operating frequency
- xMATRIX array TEE with 2,500 elements
- 2D, color flow, PW Doppler, CW Doppler, M-mode, harmonic imaging, and advanced XRES
- Electrocautery suppression
- Electronically rotating image plane from 0 to 180 degrees
- Adult transesophageal applications
- Optional adapter for use on iE33 systems

### Non-imaging









#### D2cwc CW transducer (Pedoff)

- Dedicated 2 MHz continuous wave Doppler
- Adult cardiac applications

#### D5cwc CW transducer (Pedoff)

- Dedicated 5 MHz continuous wave Doppler
- Deep venous and arterial applications

### 5.3 Transducer application guide

Transducers									
Name		C5-1	C8-5	C9-3v	L12-3	S5-1	S8-3	S12-4	X7-2t
Type of array		Curved	Curved	Curved	Linear	Sector	Sector	Sector	Sector
Number of elements		160	128	128	160	80	96	96	2500
Scanplane aperture			22.4 mm		38 mm	20.3 mm	15.4 mm	9.8 mm	Proprietary
Field of view		70°	92°	135°			90°	90°	
Broadband frequency range		5 to 1 MHz	8 to 5 MHz	9 to 3 MHz	12 to 3 MHz	5 to 1 MHz	8 to 3 MHz	12 to 4 MHz	7 to 2 MHz
Frequency									
PureWave crystal technology		●		●		●			●
Application	Exam type								
Abdominal	General	●				●			
	Renal	●				●			
	Vascular	●				●			
	Difficult	●							
Cardiac	Adult transthoracic					●	●	●	
	Adult TEE								●
	Congenital					●			
	Pediatric						●	●	
	Perioperative								●
	Intraoperative								●
Vascular	Carotid		●		●				
	Arterial		●		●				
	Venous		●		●				
	Abdominal	●				●			
	TCD					●			
	Intervention				●				
	Vascular access				●				
Obstetrics	Early OB	●							
	General OB	●		●					
	Fetal echo	●							
Gynecology	Pelvis	●		●					
	Fertility	●		●					
Small Parts	Superficial				●				
	Thyroid				●				
	Testicle				●				
	Breast				●				
Musculoskeletal	Superficial				●				
	General				●				
Acute	FAST	●				●			
	Lung	●			●				
	Vascular access				●				
Nerve	Superficial				●				
	Deep	●			●				
Pediatric	Head		●					●	
	Abdomen		●				●	●	
	Renal		●						

# 6. Measurements and analysis

## 6.1 Comprehensive measurement tools

- Cardiac analysis
- High Q automatic Doppler analysis
- QLAB quantification

## 6.2 QLAB advanced quantification

### Region of Interest (ROI) quantification plug-in

- Pixel intensity analysis, data types: echo, velocity (color) or power (angio)
- Up to 10 user-defined regions
- Thumbnail display of frames for easy trimming
- TDI velocity timing measurement
- Log/linear data display selection
- Smoothed data display option
- Vascularization Index, Flow Index and Vascularization Flow Index results on color mode files
- Motion compensation for multiframe objects

### Tissue Motion Quantification (TMQ) plug-in

- 2D cardiac quantification
- Assess global, regional/segmental, and annular (valve) cardiac function using tissue (speckle) tracking with the S5-1 transducer
  - Non-directional (Area cm<sup>2</sup>, Volume ml, Speed cm/s)
  - Directional (fractional shortening %, strain %, strain rate l/s)
  - Time motion annular displacement (TMAD) quantification

### Strain Quantification (SQ) plug-in

- For evaluation of regional myocardial function, assessment of synchronicity and guidance during bi-ventricular pacing procedures
- Tissue Doppler Imaging (TDI) velocity quantification
- Measures the myocardial velocity and derives the strain rate and strain along user-defined M-lines
- User defined M-Line motion to follow the myocardial motion
- User-selectable waveforms for optimal sub-region visualization
- Curve processing mode
- TDI velocity timing measurement

### Intima Media Thickness (IMT) measurement plug-in

- Automated assessment of the IMT on user selected frames
- For carotid and other superficial arteries

### MicroVascular Imaging (MVI) plug-in

- Integration and processing of images in contrast specific imaging mode providing detection and display of very low velocity flows of very low signal amplitude
- Motion compensation for multiframe objects

### GI 3D Quantification (GI 3DQ) plug-in

- Includes iSlice, thick slice, and MPR and volume visualization
- 3D/4D viewer for Ob/Gyn and general imaging including interventional applications
- Review 3D/4D and color 3D volumes
- Multiplanar reconstruction (MPR)
- iSlice precision volumetric slicing capability
  - Displays 2D/color slices from static or live volume
  - User-selectable slice display: 4, 9, 16 or 25
  - User-selectable interval spacing, slicing depth, and slicing source (x, y, or z)
- Free rotation of any source
- Full cine/loop review control
- 2D grayscale display adjustments
- Color display adjustments
- Zoom control
- Cine/pan slice control through volume
- User-selectable image storage
- Quick launch to measurements
  - Auto ruler display
- Compatible with freehand volumes
- Thick slice imaging
  - User-adjustable slice thickness and depth
  - Variable thick slice display adjustments with present vision settings
- 2D and 3D measurement tool including distance, area, angle, stacked and auto contour, and ellipsoid measurements
  - Vascularization Index, Flow Index and Vascularization Flow Index results on 3D color mode data sets
- Advanced XRES speckle noise reduction of MPR and volume displays
- Assisted auto-trace volume measurement tools for stacked contours and ellipse methods
- Edge detection selection for hypoechoic or high contrast targets

# 7. Physical specifications

## System dimensions

Length	14 in/35.6 cm
Width	16.25 in/41.3 cm
Height	3 in/7.6 cm
Weight	13.6 lb/6.17 kg
Display	15.4 in/39.1 cm high-resolution display with wide viewing angle

The CX50 is a laptop-sized ultrasound system you can hand carry or use on a specially designed cart.



## Physical features

- High-resolution display with wide viewing angle
- Laptop style alphanumeric QWERTY keyboard
- Ergonomic integrated carrying handle
- USB footswitch
- Multiple transducer module option multiport adaptor connects up to three transducers

## Mobility cart

- Weight: 97 lb/44 kg
- Width: 19.4 in/49.3 cm
- Depth including handle: 23.3 in/59.2 cm
- Height: adjustable from 42.5 in/107.9 cm to 32 in/82.8 cm
- Rear-mounted handle; side-mounted grips for micro-positioning
- Casters: 5 in; casters provide total locking (directional and rotational) engaged by the foot pedals
- Quick release tray
  - Simple latch system to secure CX50 system in place
  - System's integrated handle accessible from front for secure maneuverability
- Storage
  - With multiport adaptor
    - Two rear mounted storage trays
    - Cable catch tray below transducer connector
  - Without multiport adaptor ---
    - Three shelf options – two for printers, one storage, 9.9 in/25.2 cm x 13.4 in/34 cm (L x W)
    - One sliding drawer, 9.8 in/24.9 x 9.2 in/23.5 cm x 2.4 in/6 cm (L x W x D)
  - Integrated four transducer connector holders
  - Cable management hooks
- Transducer holders accommodate six transducers
- Input power to B/W printer and color printer
- USB hub with cables to:
  - B/W printer
  - Color printer
  - CX50 USB

## Travel case

- Dimensions: 22.5 in/57 cm x 16.5 in/42 cm x 12.5 in/32 cm (H x W x D)
- Weight: 15.4 lbs/7 kg when empty



- Features:

- Accommodates CX50 system and AC adapter
- Removable transducer bag stores three transducers and one gel bottle
- Wheels and retractable handle allow easy transport

#### Power requirements

System/AC adapter	100-240V, 50/60 Hz, 250 VA
System with cart and peripherals	100-240V, 50/60 Hz, 500 VA

#### Power management

- Internal lithium-ion polymer battery
- Fully charged new battery yields approx. 45 minute battery life under continuous use without AC. Actual time varies with age and condition of battery.
- Quick-charge battery technology
- Advanced battery/AC monitoring circuitry includes on-screen graphics, and low battery warning
- Suspend mode for instantaneous boot-up between exams



The easy portability of the CX50 allows premium imaging capabilities at any location, including local screening venues.

#### Environmental

Heat dissipation	700-1100 BTUs/hour (fully loaded)
Operation range	10°C–40°C operating in 15–95% relative humidity

#### ECG and physio

- One three-lead ECG input
- One external ECG input
- Two physio input channels (1V, p-p)
- Selectable ECG triggered skipping between 1 and 20
- Respiration

#### Electrical safety standards

- CSA C22.2 No. 601.1
- IEC 60601-1
- UL 60601-1
- EN 60601-1

#### Safety requirements

- Electromechanical standards met
  - C22.2 No. 601.1, Canadian Standards Association, Standard for Medical Electrical Equipment
  - JIS T 0601-1, Japanese Standard for Medical Electrical Equipment
  - EN 60601-1, European Norm, Safety of Medical Electrical Equipment
  - EN 60601-1-2 European Norm, Collateral Standard: Electromagnetic compatibility
  - EN 60601-2-37 European Norm, Particular requirements for the safety of ultrasonic medical diagnostic and monitoring equipment
  - UL 60601-1 Underwriters Laboratories Standard for Medical Electrical Equipment
- Agency approvals
  - Canadian Standards Association (CSA)
  - CE Mark in accordance with the European Medical Device Directive issued by British Standards Institute (BSI)
  - Japanese Ministry of Health, Labor and Welfare

#### Maintenance and serviceability

- Remote access for clinical and technical support\*
- Flexible service agreements
- Clinical application and educational support
- Scheduled preventative maintenance and system optimization

\* Remote Services Network – requires Philips service contract and internet access

**Philips Healthcare is part of  
Royal Philips Electronics**

How to reach us

[www.philips.com/healthcare](http://www.philips.com/healthcare)  
[healthcare@philips.com](mailto:healthcare@philips.com)

Asia

+49 7031 463 2254

Europe, Middle East, Africa

+49 7031 463 2254

Latin America

+55 11 2125 0744

North America

+1 425 487 7000

800 285 5585 (toll free, US only)

For more information visit [www.philips.com/CX50](http://www.philips.com/CX50)



© 2011 Koninklijke Philips Electronics N.V.  
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Printed in The Netherlands.  
4522 962 74831 \* JUL 2011